

Human TIM-3 protein, human IgG1 Fc tag

Cat. No. GTX65718-pro**Applications** Functional Assay**Species** Human**Package**

100 µg

Applications

Application Note

Measured by its binding ability in a functional ELISA.

Properties

Form Lyophilized powder**Buffer** Reconstitute with 100 µl ddH₂O. Lyophilized from 0.2µm-filtered PBS.**Preservative** No preservatives**Storage** Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.**Region/Sequence** The extracellular domain of human Tim-3 (aa 22-200) is fused to the N-terminus of the Fc region of human IgG1.**Expression System** CHO cells**Purity** ≥98% by SDS-PAGE.**Endotoxin** < 0.06 EU/µg**Conjugation** Unconjugated**Note** For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.For full product information, images and publications, please visit our [website](#).