

Human NCAM protein, mouse IgG2a Fc tag

Cat. No. GTX65733-pro
Applications Functional Assay

Species Human

Package

50 µg

Applications

Application Note

Measured by its binding ability in a functional ELISA.

Properties

Form Lyophilized powder

Buffer Reconstitute with 50 µl ddH₂O. Lyophilized from 0.2µm-filtered PBS.

Preservative No preservatives

Storage Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.

Region/Sequence The extracellular domain of human CD56 (a.a. 20-718) is fused to the N-terminus of the Fc region of mouse IgG2a.

Expression System HEK293 cells

Purity ??98% by SDS-PAGE.

Endotoxin < 5 EU/mg

Conjugation Unconjugated

Note For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


For full product information, images and publications, please visit our [website](#).