

Human ICOS protein, mouse IgG2a Fc tag (active)

Cat. No. GTX65736-pro
Applications Functional Assay

Species Human

Package

100 µg

Applications

Application Note

Measured by its ability to inhibit human T cell proliferation induced by Human ICOS Ligand + Human IgG1 Fc fusion protein (GTX65695-pro) in the presence of anti-CD3.

Properties

Form Lyophilized powder

Buffer Reconstitute with 100 µl ddH₂O. Lyophilized from 0.2µm-filtered PBS.

Preservative No preservatives

Storage Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.

Region/Sequence The extracellular domain of human CD278 [ICOS] (a.a. 21-139) is fused to the N-terminus of the Fc region of mouse IgG2a.

Expression System CHO cells

Purity ??98% by SDS-PAGE.

Endotoxin < 0.06 EU/µg

Conjugation Unconjugated

Note For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


For full product information, images and publications, please visit our [website](#).