

Human CD134 protein, His tag

Cat. No. GTX65756-pro

Species Human

Package

25 µg

Properties

Form Lyophilized powder

Buffer Reconstitute with ddH₂O but not less than 100 µg/ml. Lyophilized from 50mM Tris-HCl buffer (pH8.0).

Preservative No preservatives

Storage Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.

Region/Sequence C-terminal His-Tag; a.a. 29-216

Expression System E. coli

Purity ≥97% by SDS-PAGE.

Endotoxin < 0.1 EU/µg

Conjugation Unconjugated

Note For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).