

Human CD137 protein, His tag

Cat. No. GTX65757-pro

Species	Human
----------------	-------

Package

25 µg

Properties

Form	Lyophilized powder
-------------	--------------------

Buffer	Reconstitute with ddH ₂ O but not less than 100 µg/ml. Lyophilized from 50mM Tris-HCl buffer (pH8.0).
---------------	--

Preservative	No preservatives
---------------------	------------------

Storage	Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.
----------------	---

Region/Sequence	C-terminal His-Tag; a.a. 19-184
------------------------	---------------------------------

Expression System	E. coli
--------------------------	---------

Purity	≥95% by SDS-PAGE.
---------------	-------------------

Endotoxin	< 0.1 EU/µg
------------------	-------------

Conjugation	Unconjugated
--------------------	--------------

Note	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
-------------	---



For full product information, images and publications, please visit our [website](#).