

## Human CD27 protein, His tag

**Cat. No. GTX65758-pro****Species** Human**Package**

50 µg

## Properties

**Form** Lyophilized powder**Buffer** Reconstitute with ddH<sub>2</sub>O but not less than 100 µg/ml. Lyophilized from 50mM Tris-HCl buffer (pH8.0), 500mM NaCl.**Preservative** No preservatives**Storage** Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.**Region/Sequence** C-terminal His-Tag; a.a. 21-192**Expression System** E. coli**Purity** ≥97% by SDS-PAGE.**Endotoxin** < 0.1 EU/µg**Conjugation** Unconjugated**Note** For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

For full product information, images and publications, please visit our [website](#).