

Mouse CD276 protein, mouse IgG2a Fc tag

Cat. No. GTX65762-pro

Species	Mouse
----------------	-------

Package	100 µg
----------------	--------

Properties

Form	Lyophilized powder
Buffer	Reconstitute with 0.2µm-filtered PBS to 100 µg/ml. Lyophilized from 0.2µm-filtered PBS.
Preservative	No preservatives
Storage	Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.
Region/Sequence	The extracellular domain of mouse CD276 [B7-H3] (a.a. 30-248) is fused to the N-terminus of the Fc region of mouse IgG2a.
Expression System	CHO cells
Purity	≥98% by SDS-PAGE.
Endotoxin	< 0.06 EU/µg
Conjugation	Unconjugated
Note	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).