

## Mouse B7-H4 protein, mouse IgG2a Fc tag (active)

**Cat. No. GTX65763-pro**

<b>Applications</b>	Functional Assay
<b>Species</b>	Mouse

**Package**  
100 µg

## Applications

**Application Note**

Measured by its ability to inhibit anti-CD3-induced proliferation of stimulated mouse T cells.

## Properties

<b>Form</b>	Lyophilized powder
<b>Buffer</b>	Reconstitute with 100 µl ddH <sub>2</sub> O. Lyophilized from 0.2µm-filtered PBS.
<b>Preservative</b>	No preservatives
<b>Storage</b>	Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.
<b>Region/Sequence</b>	The extracellular domain of mouse B7-H4 (a.a. 29-258) is fused to the N-terminus of the Fc region of mouse IgG2a.
<b>Expression System</b>	CHO cells
<b>Purity</b>	≥98% by SDS-PAGE.
<b>Endotoxin</b>	< 0.06 EU/µg
<b>Conjugation</b>	Unconjugated
<b>Note</b>	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).