

Mouse PD-L2 protein, mouse IgG2a Fc tag

Cat. No. GTX65773-pro**Species** Mouse**Package**
100 µg

Properties

Form Lyophilized powder**Buffer** Reconstitute with 0.2µm-filtered PBS to 100 µg/ml. Lyophilized from 0.2µm-filtered PBS.**Preservative** No preservatives**Storage** Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.**Region/Sequence** The extracellular domain of mouse CD273 [PD-L2] (a.a. 20-219) is fused to the N-terminus of the Fc region of mouse IgG2a.**Expression System** CHO cells**Purity** ≥98% by SDS-PAGE.**Endotoxin** < 0.06 EU/µg**Conjugation** Unconjugated**Note** For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.For full product information, images and publications, please visit our [website](#).