

Mouse TIM-1 protein, mouse IgG2a Fc tag

Cat. No. GTX65774-pro**Species** Mouse**Package**
100 µg

Properties

Form	Lyophilized powder
Buffer	Reconstitute with 100 µl ddH ₂ O. Lyophilized from 0.2µm-filtered PBS.
Preservative	No preservatives
Storage	Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.
Region/Sequence	The extracellular domain of mouse Tim-1 (a.a. 22-212) is fused to the N-terminus of the Fc region of mouse IgG2a.
Expression System	HEK293 cells
Purity	≥98% by SDS-PAGE.
Endotoxin	< 5 EU/mg
Conjugation	Unconjugated
Note	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 10 Page 1 of 1