

Mouse TIM-3 protein, mouse IgG2a Fc tag

Cat. No. GTX65776-pro

Applications Functional Assay

Species Mouse

Package

100 µg

Applications

Application Note

Immobilized recombinant human Galectin-9 at 500ng/ml can bind mouse Tim-3 fc fusion protein with an apparent KD: <10nM. Optimal dilutions should be determined by each laboratory for each application.

Properties

Form Lyophilized powder

Buffer Reconstitute with 100 µl ddH₂O. Lyophilized from 0.2µm-filtered PBS.

Preservative No preservatives

Storage Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.

Region/Sequence The extracellular domain of mouse Tim-3 (a.a. 21-189) is fused to the N-terminus of the Fc region of mouse IgG2a.

Expression System HEK293 cells

Purity ??98% by SDS-PAGE.

Endotoxin < 5 EU/mg

Conjugation Unconjugated

Note For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).