

Mouse TIM-4 protein, mouse IgG2a Fc tag (active)

Cat. No. GTX65777-pro

Applications Functional Assay**Species** Mouse**Package**

100 µg

Applications

Application Note

Measured by its ability to inhibit anti-CD3-induced proliferation of stimulated human T cells.

Properties

Form Lyophilized powder**Buffer** Reconstitute with 100 µl ddH₂O. Lyophilized from 0.2µm-filtered PBS.**Preservative** No preservatives**Storage** Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.**Region/Sequence** The extracellular domain of mouse Tim-4 (a.a. 22-279) is fused to the N-terminus of the Fc region of mouse IgG2a.**Expression System** HEK293 cells**Purity** ??98% by SDS-PAGE.**Endotoxin** < 5 EU/mg**Conjugation** Unconjugated**Note** For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

For full product information, images and publications, please visit our [website](#).