

Mouse CD134 protein, human IgG1 Fc tag

Cat. No. GTX65778-pro

Applications Functional Assay

Species Mouse

Package

100 µg

Applications

Application Note

Measured by its binding ability in a functional ELISA.

Properties

Form Lyophilized powder

Buffer Reconstitute with 100µl ddH₂O to 1mg/ml. Lyophilized from 0.2µm-filtered PBS.

Preservative No preservatives

Storage Store at -20°C or below. After reconstitution, keep as concentrated solution. Aliquot and avoid freeze-thaw cycles.

Region/Sequence The extracellular domain of mouse CD134 (a.a. 20-211) is fused to the N-terminus of the Fc region of human IgG1.

Expression System HEK293 cells

Purity ??98% by SDS-PAGE.

Endotoxin < 5 EU/mg

Conjugation Unconjugated

Note For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).