

RanGAP1 antibody

Cat. No. GTX65813

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500 - 1:2000
ICC/IF	1:50 - 1:200
IHC-P	1:50 - 1:200

Not tested in other applications.

Calculated MW 64 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 398-587 of human RANGAP1 (NP_002874.1).
Purification	Purified by affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.




For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX65813 IHC-P Image

IHC-P analysis of human esophageal cancer tissue using GTX65813 RanGAP1 antibody.
Dilution : 1:50




GTX65813 ICC/IF Image

ICC/IF analysis of MCF7 cells using GTX65813 RanGAP1 antibody.



GTX65813 WB Image

WB analysis of various sample lysates using GTX65813 RanGAP1 antibody.
Dilution : 1:1000
Loading : 25µg per lane



GTX65813 IHC-P Image

IHC-P analysis of human placenta tissue using GTX65813 RanGAP1 antibody.
Dilution : 1:50



For full product information, images and publications, please visit our [website](https://www.genetex.com).