

SSA1 antibody

Cat. No. GTX65834

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P, IP
Reactivity	Human, Mouse, Rat

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500 - 1:2000
ICC/IF	1:50 - 1:200
IHC-P	1:50 - 1:200
IP	Assay dependent

Not tested in other applications.

Calculated MW 54 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 80-330 of human TRIM21 (NP_003132.2).
Purification	Purified by affinity chromatography
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES




GTX65834 ICC/IF Image

ICC/IF analysis of THP-1 cells using GTX65834 SSA1 antibody.

Red : Primary antibody

Blue : DAPI

Dilution : 1:200




GTX65834 IHC-P Image

IHC-P analysis of mouse spinal cord tissue using GTX65834 SSA1 antibody.

Antigen retrieval : High pressure antigen retrieval with 10 mM citrate buffer pH 6.0

Dilution : 1:50



GTX65834 IP Image

IP analysis of Jurkat cell lysates using GTX65834 SSA1 antibody.

Lane 1 : Input

Lane 2 : Ctrl IgG

Lane 3 : Primary antibody

IP antibody : 3 µg of GTX65834

Dilution : 1:1000



For full product information, images and publications, please visit our [website](https://www.genetex.com).