

HOPX antibody

Cat. No. GTX66394

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Rat

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:200 - 1:2000
IHC-P	Assay dependent

Not tested in other applications.

Calculated MW 8 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-73 of human HOPX (NP_631957.1).
Purification	Purified by affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX66394 IHC-P Image


IHC-P analysis of human skin tissue using GTX66394 HOPX antibody.

Orange : Primary antibody

Blue : DAPI

Antigen retrieval : high pressure antigen retrieval with 0.01 M citrate buffer (pH 6.0)

Dilution : 1:200



GTX66394 WB Image

WB analysis of rat heart tissue lysate using GTX66394 HOPX antibody.

Dilution : 1:600

Loading : 25µg per lane



For full product information, images and publications, please visit our [website](#).