

PEA15 (phospho Ser104) antibody

Cat. No. GTX66602

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Human, Rat

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500 - 1:2000

Not tested in other applications.

Calculated MW 15 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	A synthetic phosphorylated peptide around S104 of human PEA15 (NP_003759.1).
Purification	Purified by affinity chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES




GTX66602 WB Image

WB analysis of HeLa cells treated by PMA/TPA (200nM) for 15 minutes after serum-starvation overnight using GTX66602 PEA15 (phospho Ser104) antibody.

Dilution : 1:2000

Loading : 25µg per lane



GTX66602 WB Image

WB analysis of HeLa cells treated by PMA/TPA (200nM) for 15 minutes after serum-starvation overnight and C6 cell lysate using GTX66602 PEA15 (phospho Ser104) antibody.

Dilution : 1:2000

Loading : 25µg per lane



For full product information, images and publications, please visit our [website](https://www.genetex.com).