

PTPRO antibody

Cat. No. GTX66615

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	IHC-P
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
IHC-P	1:50-1:200


Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant full length Human PTPRO.
Purification	Purified by antigen-affinity chromatography From serum
Purity	> 95% (by SDS-PAGE)
Conjugation	Unconjugated
Note	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES

GTX66615 IHC-P Image

IHC-P analysis of human stomach carcinoma tissue using GTX66615 PTPRO antibody.

Left : Primary antibody

Right : Negative control

Dilution : 1:50


GTX66615 IHC-P Image

IHC-P analysis of human kidney carcinoma tissue using GTX66615 PTPRO antibody.

Left : Primary antibody

Right : Negative control

Dilution : 1:50



For full product information, images and publications, please visit our [website](#).