

EphB4 antibody

Cat. No. GTX66654

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse, Rat

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:1000
IHC-P	1:50-1:200

Not tested in other applications.

Calculated MW 108 kDa. ([Note](#))


Properties

Form	Liquid
Buffer	PBS
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide, corresponding to amino acids 580-630 of Human EphB4.
Purification	Purified by antigen-affinity chromatography From serum
Purity	> 95% (by SDS-PAGE)
Conjugation	Unconjugated
Note	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX66654 WB Image

WB analysis of various samples using GTX66654 EphB4 antibody.


Lane1 : Hcc827 whole cell lysate

Lane2 : LO2 whole cell lysate

Lane3 : The Kidney tissue lysate of Mouse

Lane4 : The Kidney tissue lysate of Rat

Dilution : 1:500



GTX66654 IHC-P Image

IHC-P analysis of human brain tissue using GTX66654 EphB4 antibody.



For full product information, images and publications, please visit our [website](https://www.genetex.com).