

CDC20 antibody

Cat. No. GTX66665

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB
Reactivity	Human, Mouse, Rat

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:1000

Not tested in other applications.


Calculated MW 55 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide, corresponding to amino acids 80-130 of Human Cdc20.
Purification	Purified by antigen-affinity chromatography From serum
Purity	> 95% (by SDS-PAGE)
Conjugation	Unconjugated
Note	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX66665 WB Image

WB analysis of various samples using GTX66665 CDC20 antibody.


Lane1 : HEK293T whole cell lysate(40ug)

Lane2 : PC3 whole cell lysate(40ug)

Lane3 : The Brain tissue lysate of Mouse(40ug)

Lane4 : The Prostate tissue lysate of Rat(40ug)

Dilution : 1:500


GTX66665 WB Image


WB analysis of various samples using GTX66665 CDC20 antibody.

Lane1 : HeLa cell lysate

Lane2 : Raw264.7 cell lysate

Lane3 : Rat liver tissue lysate

Dilution : 1:500


GTX66665 WB Image

WB analysis of various samples using GTX66665 CDC20 antibody.

Lane1 : SGC7901 whole cell lysate(40ug)

Lane2 : HeLa whole cell lysate(40ug)

Lane3 : The Heart tissue lysate of Mouse(40ug)

Lane4 : The Heart tissue lysate of Rat(40ug)

Dilution : 1:500



For full product information, images and publications, please visit our [website](http://www.genetex.com).