

KSR2 antibody

Cat. No. GTX66673

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse, Rat

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:1000
IHC-P	1:50-1:200

Not tested in other applications.

Calculated MW 108 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide, corresponding to amino acids 662-724 of Human Ksr2.
Purification	Purified by antigen-affinity chromatography From serum
Purity	> 95% (by SDS-PAGE)
Conjugation	Unconjugated
Note	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX66673 IHC-P Image

IHC-P analysis of human brain tissue using GTX66673 KSR2 antibody.



GTX66673 WB Image

WB analysis of various samples using GTX66673 KSR2 antibody.


Lane1 : HEK293T whole cell lysate(40μg)

Lane2 : U-87MG whole cell lysate(40μg)

Lane3 : C6 whole cell lysate(40μg)

Lane4 : The Brain tissue lysate of Mouse (40μg)

Dilution : 1:1000



GTX66673 WB Image

WB analysis of various samples using GTX66673 KSR2 antibody.

Lane1 : HEK293T whole cell lysate(40ug)

Lane2 : A549 whole cell lysate(40ug)

Lane3 : AML-12 whole cell lysate(40ug)

Lane4 : C6 whole cell lysate(40ug)

Dilution : 1:500



For full product information, images and publications, please visit our [website](https://www.genetex.com).