

SFRS3 antibody

Cat. No. GTX66685

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse, Rat

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:1000
IHC-P	1:50-1:200


Not tested in other applications.

Calculated MW 19 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide, corresponding to amino acids 120-170 of Human SRp20.
Purification	Purified by antigen-affinity chromatography From serum
Purity	> 95% (by SDS-PAGE)
Conjugation	Unconjugated
Note	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX66685 IHC-P Image

IHC-P analysis of human liver carcinoma tissue using GTX66685 SFRS3 antibody.

Left : Primary antibody

Right : Negative control

Dilution : 1:50


GTX66685 WB Image

WB analysis of various samples using GTX66685 SFRS3 antibody.


Lane1 : U-87MG whole cell lysate(40ug)

Lane2 : Jurkat whole cell lysate(40ug)

Lane3 : C6 whole cell lysate(40ug)

Lane4 : BV2 whole cell lysate(40ug)

Dilution : 1:500


GTX66685 IHC-P Image

IHC-P analysis of human stomach carcinoma tissue using GTX66685 SFRS3 antibody.

Left : Primary antibody

Right : Negative control

Dilution : 1:50



For full product information, images and publications, please visit our [website](https://www.genetex.com).