

IRF4 / MUM1 antibody

Cat. No. GTX66692

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse, Rat

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-P	1:50-1:200

Not tested in other applications.


Calculated MW 52 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	15mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide, corresponding to amino acids 262-318 of Human IRF4 / MUM1.
Purification	Purified by antigen-affinity chromatography From serum
Purity	> 95% (by SDS-PAGE)
Conjugation	Unconjugated
Note	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX66692 WB Image


WB analysis of various samples using GTX66692 IRF4 / MUM1 antibody.

Lane1 : HEK293T whole cell lysate

Lane2 : Raw264.7 whole cell lysate

Lane3 : PC12 whole cell lysate

Dilution : 1:500



GTX66692 IHC-P Image

IHC-P analysis of human tonsil carcinoma tissue using GTX66692 IRF4 / MUM1 antibody.

Left : Primary antibody

Right : Negative control

Dilution : 1:50



For full product information, images and publications, please visit our [website](https://www.genetex.com).