

PMS2 antibody

Cat. No. GTX66694

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse, Rat

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
IHC-P	1:50-1:200

Not tested in other applications.


Calculated MW 96 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide, corresponding to amino acids 450-500 of Human PMS2.
Purification	Purified by antigen-affinity chromatography From serum
Purity	> 95% (by SDS-PAGE)
Conjugation	Unconjugated
Note	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).


DATA IMAGES

GTX66694 IHC-P Image

IHC-P analysis of human kidney carcinoma tissue using GTX66694 PMS2 antibody.

Left : Primary antibody

Right : Negative control

Dilution : 1:50


GTX66694 WB Image

WB analysis of various samples using GTX66694 PMS2 antibody.


Lane1 : A549 whole cell lysate

Lane2 : Raw264.7 whole cell lysate

Lane3 : H9C2 whole cell lysate

Lane4 : HeLa whole cell lysate

Dilution : 1:500



GTX66694 IHC-P Image

IHC-P analysis of human colon carcinoma tissue using GTX66694 PMS2 antibody.

Left : Primary antibody

Right : Negative control

Dilution : 1:50


GTX66694 WB Image

WB analysis of various samples using GTX66694 PMS2 antibody.

Lane1 : A549 whole cell lysate

Lane2 : Raw264.7 whole cell lysate

Lane3 : H9C2 whole cell lysate

Dilution : 1:500



For full product information, images and publications, please visit our [website](https://www.genetex.com).