

Interferon gamma antibody

Cat. No. GTX66714

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse, Rat

References (4)

Package

100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:1000
IHC-P	1:50-1:200


Not tested in other applications.

Calculated MW 19 kDa. ([Note](#))


Properties

Form	Liquid
Buffer	PBS
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide, corresponding to amino acids 31-80 of Human IFN-γ.
Purification	Purified by antigen-affinity chromatography From serum
Purity	> 95% (by SDS-PAGE)
Conjugation	Unconjugated
Note	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX66714 IHC-P Image

IHC-P analysis of human lymph node tissue using GTX66714 Interferon gamma antibody.


GTX66714 WB Image


WB analysis of various samples using GTX66714 Interferon gamma antibody.

Lane1 : HEK293T whole cell lysate

Lane2 : Raw264.7 whole cell lysate

Lane3 : H9C2 whole cell lysate

Dilution : 1:500


GTX66714 WB Image

WB analysis of various samples using GTX66714 Interferon gamma antibody.

Lane1 : Jurkat whole cell lysate(40ug)

Lane2 : A549 whole cell lysate(40ug)

Lane3 : The spleen tissue lysate of Mouse(30ug)

Lane4 : The spleen tissue lysate of Rat(30ug)

Lane5 : SP2/0 whole cell lysate

Dilution : 1:500



For full product information, images and publications, please visit our [website](https://www.genetex.com).