

PSMD11 antibody

Cat. No. GTX66726

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse, Rat

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:1000
IHC-P	1:50-1:200

Not tested in other applications.

Calculated MW 47 kDa. ([Note](#))

Properties


Form	Liquid
Buffer	PBS
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic peptide, corresponding to amino acids 281-330 of Human PSMD11.
Purification	Purified by antigen-affinity chromatography From serum
Purity	> 95% (by SDS-PAGE)
Conjugation	Unconjugated
Note	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 06 Page 1 of 2

DATA IMAGES




GTX66726 IHC-P Image

IHC-P analysis of human breast carcinoma tissue using GTX66726 PSMD11 antibody.
Dilution : 1:100




GTX66726 WB Image

WB analysis of various samples using GTX66726 PSMD11 antibody.
Lane1 : HCC827 whole cell lysate(40ug)
Lane2 : HCT116 whole cell lysate(40ug)
Lane3 : 786-O whole cell lysate(40ug)
Lane4 : The Muscle tissue lysate of Mouse(40ug)
Lane5 : PC12 whole cell lysate(40ug)
Dilution : 1:500



GTX66726 WB Image

WB analysis of various samples using GTX66726 PSMD11 antibody.
Lane1 : HEK293T whole cell lysate(40ug)
Lane2 : PANC1 whole cell lysate(40ug)
Lane3 : L02 whole cell lysate(40ug)
Lane4 : H9C2 whole cell lysate(40ug)
Lane5 : CT26 whole cell lysate(40ug)
Dilution : 1:500



GTX66726 WB Image

WB analysis of various samples using GTX66726 PSMD11 antibody.
Lane1 : A549 whole cell lysate(40ug)
Lane2 : K562 whole cell lysate(40ug)
Lane3 : The Heart tissue lysate of Rat(40ug)
Lane4 : The Heart tissue lysate of Mouse(40ug)
Dilution : 1:500



For full product information, images and publications, please visit our [website](https://www.genetex.com).