

PKC zeta (phospho Thr410) antibody

Cat. No. GTX66730

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, IHC-P
Reactivity	Human, Mouse, Rat

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:1000
IHC-P	1:50-1:200


Not tested in other applications.

Calculated MW 68 kDa. ([Note](#))


Properties

Form	Liquid
Buffer	PBS
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic phosphopeptide derived from human PKC ζ around the phosphorylation site of Threonine 410.
Purification	Purified by antigen-affinity chromatography From serum
Purity	> 95% (by SDS-PAGE)
Conjugation	Unconjugated
Note	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX66730 IHC-P Image

IHC-P analysis of human breast carcinoma tissue using GTX66730 PKC zeta (phospho Thr410) antibody.
Dilution : 1:100


GTX66730 WB Image

WB analysis of various samples using GTX66730 PKC zeta (phospho Thr410) antibody.

Lane1 : MCF-7 whole cell lysate

Lane2 : MCF-7 treated with insulin(100nM/ml,15 minutes) whole cell lysate

Lane3 : The Kidney tissue lysate of Rat

Lane4 : The Lung tissue lysate of Mouse

Dilution : 1:500



For full product information, images and publications, please visit our [website](https://www.genetex.com).