

## Adrenomedullin antibody

Cat. No. GTX66817

|              |            |
|--------------|------------|
| Host         | Rabbit     |
| Clonality    | Polyclonal |
| Isotype      | IgG        |
| Applications | WB         |
| Reactivity   | Mouse, Rat |

Package  
100 µl

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | 1:500-1:1000         |

Not tested in other applications.

Calculated MW 20 kDa. ( [Note](#) )


## Properties

|               |                                                                                                                                                                                                                            |
|---------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Form          | Liquid                                                                                                                                                                                                                     |
| Buffer        | PBS                                                                                                                                                                                                                        |
| Preservative  | 0.05% Sodium azide                                                                                                                                                                                                         |
| Storage       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| Concentration | 1 mg/ml (Please refer to the vial label for the specific concentration.)                                                                                                                                                   |
| Immunogen     | Synthetic peptide, corresponding to amino acids 116-160 of Human ADM                                                                                                                                                       |
| Purification  | Purified by antigen-affinity chromatography<br>From serum                                                                                                                                                                  |
| Purity        | > 95% (by SDS-PAGE)                                                                                                                                                                                                        |
| Conjugation   | Unconjugated                                                                                                                                                                                                               |
| Note          | For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.                                                                          |



For full product information, images and publications, please visit our [website](#).

## DATA IMAGES



### GTX66817 WB Image

WB analysis of various samples using GTX66817 Adrenomedullin antibody.

Lane1 : HCT116 whole cell lysate(40ug)

Lane2 : A549 whole cell lysate(40ug)

Lane3 : HEK293T whole cell lysate(40ug)

Lane4 : The Kidney lysate of Rat(40ug)

Lane5 : The Kidney lysate of Mouse (40ug)

Dilution : 1:500



For full product information, images and publications, please visit our [website](https://www.genetex.com).