

BRG1 antibody [6D7]

Cat. No. GTX66842

Host	Mouse
Clonality	Monoclonal
Isotype	IgG
Applications	WB, ICC/IF, IP
Reactivity	Human, Mouse

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000
ICC/IF	1:50-1:200
IP	1:50-1:200

Not tested in other applications.

Calculated MW 185 kDa. ([Note](#))

Properties


Form	Liquid
Buffer	100mM Tris-Glycine, 150mM NaCl, 50% Glycerol
Preservative	0.2% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Purified recombinant human SMARCA4 protein fragments expressed in E.coli.
Purification	Purified by antigen-affinity chromatography
Purity	> 95% (by SDS-PAGE)
Conjugation	Unconjugated
Note	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).

Date 2026 / 10 / 07 Page 1 of 2

DATA IMAGES




GTX66842 ICC/IF Image

ICC/IF analysis of 4% PFA-fixed HeLa cells using GTX66842 BRG1 antibody [6D7].


Blue : DAPI

Dilution : 1:50



GTX66842 IP Image

IP analysis of K562 cell lysates using GTX66842 BRG1 antibody [6D7].



GTX66842 WB Image

WB analysis of HeLa NE, F9, and K562 cell lysates using GTX66842 BRG1 antibody [6D7].

Dilution : 1:1000



For full product information, images and publications, please visit our [website](https://www.genetex.com).