

RAD23B antibody [5H1]

Cat. No. GTX66845

Host	Mouse
Clonality	Monoclonal
Isotype	IgG
Applications	WB, ICC/IF, IHC-P
Reactivity	Human, Mouse, Hamster, Monkey

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:1000
ICC/IF	1:50-1:200
IHC-P	1:50-1:200

Not tested in other applications.

Calculated MW 43 kDa. ([Note](#))

Properties


Form	Liquid
Buffer	PBS, 50% Glycerol
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Purified recombinant human Rad23B protein fragments expressed in E.coli.
Purification	Purified by antigen-affinity chromatography
Purity	> 95% (by SDS-PAGE)
Conjugation	Unconjugated
Note	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).


Date 2026 / 10 / 07 Page 1 of 2

DATA IMAGES




GTX66845 WB Image

WB analysis of K562, SW480, CHO-K1, 3T3, and COS7 cell lysates using GTX66845 RAD23B antibody [5H1].
Dilution : 1:1000



GTX66845 IHC-P Image

IHC-P analysis of prostate cancer tissue using GTX66845 RAD23B antibody [5H1].
Antigen retrieval : citrate buffer (pH 6.0)



GTX66845 ICC/IF Image

ICC/IF analysis of 4% PFA-fixed HeLa cells using GTX66845 RAD23B antibody [5H1].
Dilution : 1:100



For full product information, images and publications, please visit our [website](https://www.genetex.com).