

CD11b antibody [L1]

Cat. No. GTX66849

Host	Mouse
Clonality	Monoclonal
Isotype	IgG
Applications	WB
Reactivity	Human, Mouse, Rat

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:1000

Not tested in other applications.


Calculated MW 127 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	100mM Tris-Glycine, 150mM NaCl, 50% Glycerol
Preservative	0.2% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant Human CD11B (630aa-1104aa).
Purification	Purified by antigen-affinity chromatography
Purity	> 95% (by SDS-PAGE)
Conjugation	Unconjugated
Note	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX66849 WB Image

WB analysis of various samples using GTX66849 CD11b antibody [L1].

Lane1 : Myla2059 whole cell lysate(30ug)

Lane2 : HuT78 whole cell lysate(30ug)

Lane3 : K562 whole cell lysate(30ug)

Lane4 : The Spleen tissue lysate of Rat(40ug)

Lane5 : RAW264.7 whole cell lysate(40ug)

Lane6 : THP-1 whole cell lysate(40ug)

Dilution : 1:500



For full product information, images and publications, please visit our [website](https://www.genetex.com).