

Human ACP1 protein, His tag (active)

Cat. No. GTX66900-pro

Applications Functional Assay**Species** Human**Package**

100 µg

Applications

Application Note


Specific activity is > 15000 units/mg, and is defined as the amount of enzyme that hydrolyze 1.0 nmoles of p-nitrophenyl phosphate (pNPP) per minute at pH 7.5 at 37°C.

Properties

Form Liquid**Buffer** 20mM MES buffer, 10% Glycerol, 0.1mM PMSF, 2mM EDTA**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 1 mg/ml (Please refer to the vial label for the specific concentration.)**Region/Sequence** Full length protein, N-terminal His-Tag; Length: 178 a.a. Sequence: MGSSHHHHHH SSGLVPRGSH MAEQATKSVL FVCLGNICRS PIAEAVFRKL VTDQNISENW VIDSGAVSDW NVGRSPDPRA VSCLRNHGIH TAHKARQITK EDFATFDYIL CMDESNLRLD NRKSNQVKTC KAKIELLSY DPQKQLIED PYYGNDSDFE TVYQQCVRCC RAFLEKAH**Expression System** E. coli**Purity** > 95% by SDS-PAGE.**Endotoxin** < 1.0 EU/µg (determined by LAL method)**Conjugation** Unconjugated**Note** For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX66900-pro Image

3 µg of GTX66900-pro Human ACP1 protein (active) by SDS-PAGE under reducing condition and visualized by coomassie blue stain



For full product information, images and publications, please visit our [website](https://www.genetex.com).