

Human Adenylate kinase 1 protein, His tag (active)

Cat. No. GTX66905-pro**Applications** Functional Assay**Species** Human**Package**

100 µg

Applications

Application Note

Specific activity: > 600 units/mg. One unit will convert 2.0 µmoles of ADP to ATP + AMP per minute at pH 7.5 at 37°C.

Properties


Form Liquid**Buffer** 20mM Tris-HCl, 10% Glycerol**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 1 mg/ml (Please refer to the vial label for the specific concentration.)**Region/Sequence**

Full length protein, N-terminal His-Tag; Length: 214 a.a. Sequence: MGSSHHHHHH SSGLVPRGSH MEEKLKTKI
IFVGGPGSG KGTQCEKIVQ KYGYTHLSTG DLLRSEVSSG SARGKKLSEI MEKQLVPLE TVLDMLRDAM VAKVNTSKGF
LIDGYPREVQ QGEFERRIG QPTLLLYVDA GPETMTQRLL KRGETSGRVD DNEETIKRRL ETTYKATEPV IAFYEKRGIV RKVNAEGSVD
SVFSQVCTHL DALK

Expression System E. coli**Purity** > 95% by SDS-PAGE.**Conjugation** Unconjugated**Note** For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX66905-pro Image

3 μ g of GTX66905-pro Human Adenylate kinase 1 protein (active) by SDS-PAGE under reducing condition and visualized by coomassie blue stain



For full product information, images and publications, please visit our [website](https://www.genetex.com).