

Human Catalase protein, His tag (active)

Cat. No. GTX66938-pro

Applications Functional Assay**Species** Human**Package**

100 µg

Applications

Application Note

Specific activity is > 70000 units/mg. One unit will decompose 1.0µmole of H₂O₂ per minute at pH 8.0 at 25°C.

Properties


Form Liquid**Buffer** 20mM Tris-HCl, 150mM NaCl, 10% Glycerol**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 1 mg/ml (Please refer to the vial label for the specific concentration.)**Region/Sequence**

Full length protein, N-terminal His-Tag; Length: 547 a.a. Sequence: MGSSHHHHHH SGLVPRGSH MADSRDPASD QMQHWKEQRA AQKADVLTTG AGNPVGDKLN VITVGPRGPL LVQDVVFTDE MAHFDRERIP ERVHAKGAG AFGYFEVTHD ITKYSKAKVF EHIGKKTPIA VRFSTVAGES GSADTVRDPR GFAVKFYTED GNWDLVGNNT PIFFIRDPII FPSFIHSQKR NPQTHLKDPD MVWDFWLSLRP ESLHQVSFLF SDRGIPDGHR HMNGYGSHTF KLVNANGEAV YCKFHYKTDQ GIKNLSVEDA ARLSQEDPDY GIRDLFNAIA TGKYPSWTFY IQVMTFNQAE TFPFNPFDLT KWWPHKDYPL IPVGKLVNLR NPVNYFAEVE QIAFDPSNMP PGIEASPDKM LQGRLFAYPD THRHRLGPNY LHIPVNCYPY ARVANYQRDG PMCMQDNQGG APNYYPNSFG APEQQPSALE HSIQYSGEVR RFNTANDDNV TQVRAFVNV LNEEQRKRLC ENIAGHLKDA QIFIQKKAVK NFEVHPDYG SHIQALLDKY NAEKPKNAIH TFVQSGSHLA AREKANL

Expression System E. coli**Purity** > 90% by SDS-PAGE.**Conjugation** Unconjugated**Note** For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX66938-pro Image

3 μ g of GTX66938-pro Human Catalase protein (active) by SDS-PAGE under reducing condition and visualized by coomassie blue stain



For full product information, images and publications, please visit our [website](#).