

Human CNTF protein (active)

Cat. No. GTX66947-pro

Applications	Functional Assay
Species	Human

Package
20 µg

Applications

Application Note

Measured in a cell proliferation assay using TF-1 human erythroleukemic cells. The ED₅₀ for this effect is less than or equal to 500 ng/ml.


Properties

Form	Liquid
Buffer	20mM Tris-HCl, 10% Glycerol, 1mM DTT
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Region/Sequence	Full length protein, Length: 200 a.a. Sequence: MAFTEHSPLT PHRRDLCSRS IWLARKIRSD LTALTESYVK HQGLNKNINL DSADGMPVAS TDQWSELTEA ERLQENLQAY RTFHVLLARL LEDQQVHFTP TEGDFHQAIH TLLQVAAFA YQIEELMILL EYKIPRNEAD GMPINVGDDG LFEKKLWGLK VLQELSQWTV RSIHDLRFIS SHQTGIPARG SHYIANNKKM
Expression System	E. coli
Purity	> 90% by SDS-PAGE.
Endotoxin	< 1.0 EU/µg (determined by LAL method)
Conjugation	Unconjugated
Note	For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX66947-pro Image

3 µg of GTX66947-pro Human CNTF protein (active) by SDS-PAGE under reducing condition and visualized by coomassie blue stain



For full product information, images and publications, please visit our [website](#).