

Human IL5 protein, His tag (active)

Cat. No. GTX67068-pro

Applications Functional Assay**Species** Human**Package**

10 µg

Applications

Application Note


Measured in a cell proliferation assay using TF-1 human erythroleukemic cells. The ED₅₀ for this effects is less or equal to 1.5 ng/ml.

Observed MW (kDa) 13.5-18 kDa.

Properties

Form Liquid**Buffer** PBS, 10% Glycerol**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 0.5 mg/ml (Please refer to the vial label for the specific concentration.)**Region/Sequence** C-terminal His-Tag; Length: 124 a.a. Sequence: ADPIPTIPT SALVKETLAL LSTHRTLIA NETLRIPVPV HKNHQLCTEE IFQGIGTLES QTVQGGTVER LFKNLSLIK YIDGQKKKCG EERRRVNQFL DYLQEF LGVM NTEWIESHH HHHH**Expression System** Baculovirus**Purity** > 90% by SDS-PAGE.**Endotoxin** < 1.0 EU/µg (determined by LAL method)**Conjugation** Unconjugated**Note** For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX67068-pro Image

3 µg of GTX67068-pro Human IL5 protein (active) by SDS-PAGE under reducing condition and visualized by coomassie blue stain



For full product information, images and publications, please visit our [website](#).