

Human PGAM2 protein, His tag (active)

Cat. No. GTX67112-pro
Applications Functional Assay

Species Human

Package

20 µg

Applications

Application Note

Specific activity is > 100 units/mg, in which One unit will convert 1.0 µmole of 3-phosphoglycerate to 2-phosphoglycerate per minute at pH 7.6 at 37°C.

Properties

Form Liquid

Buffer 20mM Tris-HCl, 100mM NaCl, 20% Glycerol, 1mM DTT

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 1 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence

Full length protein, N-terminal His-Tag; Length: 273 a.a. Sequence: MGSSHHHHHH SSGLVPRGSH MATHRLVMVR HGESTWNQEN RFCGWFDAEL SEKGTEEAKR GAKAIKDAKM EFDICYTSL KRAIRTLWAI LDGTDQMWWLP VRTWRLNER HYGGLTGLNK AETAAGHGE QVKIWRRSFD IPPPPMDEKH PYNSISKER RYAGLKPGE PTCESLKDTI ARALPFWNEE IVPQIKAGKR VLIAAHGNSL RGIVKHLEGM SDQAIMELNL PTGIPIVEL NKEKPTKPM QFLGDEETVR KAMEAVAAQG KAK

Expression System E. coli


Purity > 95% by SDS-PAGE.

Conjugation Unconjugated

Note For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.


For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX67112-pro Image

3 µg of GTX67112-pro Human PGAM2 protein (active) by SDS-PAGE under reducing condition and visualized by coomassie blue stain



For full product information, images and publications, please visit our [website](https://www.genetex.com).