

Human PHPT1 protein, His tag (active)

Cat. No. GTX67120-pro

Applications Functional Assay**Species** Human**Package**

20 µg

Applications

Application Note


Specific activity is > 120 units/mg, and is defined as the amount of enzyme that hydrolyze 1nmol p-nitrophenyl phosphate per minute at pH 7.5 at 37°C.

Properties

Form Liquid**Buffer** 20mM Tris-HCl, 200mM NaCl, 10% Glycerol, 2mM DTT**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 0.5 mg/ml (Please refer to the vial label for the specific concentration.)**Region/Sequence** Full length protein, N-terminal His-Tag; Length: 145 a.a. Sequence: MGSSHHHHHH SGLVPRGSH MAVADLALIP
DVDIDSDGVF KYVLRVHSA PRSGAPAAES KEIVRGYKWA EYHADIYDKV SGDMQKQGCD CECLGGGRIS HQSQDKKIHV
YGYSMAYGPA QHAISTEKIK AKYPDYEVTW ANDGY**Expression System** E. coli**Purity** > 95% by SDS-PAGE.**Conjugation** Unconjugated**Note** For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX67120-pro Image

3 μ g of GTX67120-pro Human PHPT1 protein (active) by SDS-PAGE under reducing condition and visualized by coomassie blue stain



For full product information, images and publications, please visit our [website](#).