

Human TXN2 protein (active)

Cat. No. GTX67177-pro

Applications Functional Assay

Species Human

Package

100 µg

Applications

Application Note

Specific activity is > 40 A650/cm/min/mg, obtained by measuring the increase of insulin precipitation in absorbance at 650 nm resulting from the reduction of insulin.

Properties

Form Liquid

Buffer PBS

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 1 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence Length: 108 a.a. Sequence: MTTFNIQDGP DFQDRVNSE TPVVDFHAQ WCGPCKILGP RLEKMTAKQH GKVVMAKVDI DDHTDLAIEY EVSAVPTVLA MKNGDWDKF VGKDEDQLE AFLKKLIG

Expression System E. coli

Purity > 95% by SDS-PAGE.


Conjugation Unconjugated

Note For In vitro laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX67177-pro Image

3 μ g of GTX67177-pro Human TXN2 protein (active) by SDS-PAGE under reducing condition and visualized by coomassie blue stain



For full product information, images and publications, please visit our [website](#).