

# Human ADSL protein, His tag

**Cat. No. GTX67201-pro**
**Species** Human

**Package**  
100 µg

**Applications**
**Observed MW (kDa)** 55 kDa.

**Properties**
**Form** Liquid

**Buffer** 20mM Tris-HCl, 100mM NaCl, 40% Glycerol, 1mM DTT

**Preservative** No preservatives

**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

**Concentration** 1 mg/ml (Please refer to the vial label for the specific concentration.)

**Region/Sequence**

Full length protein, N-terminal His-Tag; MRGSHHHHHH GMASMTGGQQ MGRDLYDDDD KDRWGSMAG GDHGSPDSYR  
 SPLASRYASP EMCVFVSDRY KFRTWRQLWL WLAEAEQTLG LPITDEQIQE MKNLENIDF KMAAEEKRL RHDVMAHVHT  
 FGHCCPKAAG IHLGATSCY VGDNTDLIIL RNALDLLPK LARVISRLAD FAKERASLPT LGFTHFQPAQ LTTVGKRCCL WIQDLCMDLQ  
 NLKRVRRDLR FRGVKGTGT QASFLQLFEG DDHKVEQLDK MVTEKAGFKR AFIITGQTYT RKVDIEVLSV LASLGASVHK  
 ICTDIRLLAN LKEMEPEFEK QQIGSSAMPY KRNPMSERC CSLARHMLTL VMDPLQTASV QWFERTLDDS ANRRICLAEA  
 FLTADTILNT LQNISEGLVV YPKVIERRIR QELPFMATEN IIMAMVKAGG SRQDCHEKIR VLSQQAASVV KQEGGDNDLI  
 ERIQVDAYFS PIHSQLDHLL DPSSFTGRAS QQVQRFLEEE VYPLLPYES VMKVKAELCL

**Expression System** E. coli

**Purity** > 95% by SDS-PAGE.

**Conjugation** Unconjugated

**Note**


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

## DATA IMAGES



### GTX67201-pro Image

3μg Human ADSL protein (GTX67201-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).