

Human CDC25A protein, His tag

Cat. No. GTX67279-pro

Species Human**Package**

10 µg

Applications

Observed MW (kDa) 59 kDa.

Properties

Form Liquid**Buffer** 20mM Tris-HCl, 200mM NaCl, 20% Glycerol, 5mM DTT, 1mM EDTA**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 0.25 mg/ml (Please refer to the vial label for the specific concentration.)**Region/Sequence**

Full length protein, N-terminal His-Tag; MRGSHHHHHH GMASMTGGQQ MGRDLYDDDD KDRWGSMELG PEPHRRRL
FACSPPPASQ PVVKALFGAS AAGGLSPVTN LVTMDQLQG LGSDYEQPLE VKNNSNLQRM GSSESTDGSGF CLDSPGPLDS
KENLENPMRR IHSLPQKLLG CSPALKRSHS DSLDHDIFQL IDPDENKENE AFEFKKPVPR VSRGCLHSHG LQEGKDLFTQ
RQNSAPARML SSNERDSSEP GNFIPLFTPQ SPVTATLSDE DDGFVDLLDG ENLKNEETP SCMASLWTAP LVMRTTNLDN
RCKLFDPSPL CSSSTRSVLK RPERSQEESE PGSTKRRKSM SGASPKESTN PEKAHETLHQ SLSLASSPKG TIENILDNDP RDLIGDFSKG
YLFHTVAGKH QDLKYISPEI MASVLNGKFA NLIKEFVIID CRYPYEYEGG HIKGAVNLHM EEEVEDFLK KPIVPTDGKR VIVVFHCEFS
SERGPRMCRY VRERDRLGNE YPKLHYPELY VLKGGYKEFF MKCQSYCEPP SYRPMHHEDF KEDLKKFRTK SRTWAGEKSK
REMSRLKKL

Expression System E. coli**Purity** > 80% by SDS-PAGE.**Conjugation** Unconjugated**Note**


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX67279-pro Image

3μg Human CDC25A protein (GTX67279-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).