

Human CtBP1 protein

Cat. No. GTX67342-pro
Species Human

Package
100 µg

Applications
Observed MW (kDa) 48 kDa.

Properties
Form Liquid

Buffer 20mM Tris-HCl, 100mM NaCl, 10% Glycerol

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 1 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence

Full length protein, MGSSHLLNKG LPLGVRPPIM NGPLHPRPLV ALLDGRDCTV EMPILKDVAT VAFCD AQSTQ EIHEKVLNEA
VGALMYHTIT LTREDLEKFK ALRIIVRIGS GFDNIDIKSA GDLGIAVCNV PAASVEETAD STLCHILNLY RRATWLHQAL REGTRVQSVE
QIREVASGAA RIRGETLGII GLGRVGQAVA LRAKAFGFNV LFYDPYLSDG VERALGLQRV STLQDLLFHS DCVTLHCGLN
EHNHHLINDF TVKQMRQGAF LVNTARGGLV DEKALAQALK EGRIIGAALD VHESEPFSS QGPLKDAPNL ICTPHAAWYS
EQASIEMREE AAREIRRAIT GRIPDSLKNC VNKDHLTAAT HWASMDPAVV HPELNGAAYR YPPGVWGVAP TGIPAAVEGI
VPSAMSLSHG LPPVAHPPHA PSPGQTVKPE ADRDHASDQL

Expression System E. coli

Purity > 90% by SDS-PAGE.

Conjugation Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX67342-pro Image

3μg Human CtBP1 protein (GTX67342-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](https://www.genetex.com).