

Human CTH protein, His tag


Cat. No. GTX67343-pro

Species	Human	Package 100 µg
Applications		
Observed MW (kDa)	45 kDa.	
Properties		
Form	Liquid	
Buffer	20mM Tris-HCl, 100mM NaCl, 10% Glycerol, 2mM DTT	
Preservative	No preservatives	
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.	
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)	
Region/Sequence	Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MQEKDASSQG FLPHFQHFAT QAIHVGQDPE QWTSRAVVPP ISLSTTFKQG APGQHS GFY SRSGNPTRNC LEKAVAALDG AKYCLAFASG LAATVTITHL LKAGDQICM DDVYGGTNRY FRQVASEFGL KISFVDCSKI KLEAAITPE TKLVWETPT NPTQKVIDIE GCAHIVHKHG DIILVDNTF MSPYFQRPLA LGADISMYSY TKYMNGHSDV VMGLVSVNCE SLHNRLRFLQ NSLGAVPSPI DCYLCNRLGK TLHVRMEKHF KNGMAVAQFL ESNPWVEKVI YPGLPSHPQH ELVKRQCTGC TGMVTFYIKG TLQHAEIFLK NLKFLTAES LGGFESLAE PAIMTHASVL KNDRDVLGIS DTLIRLSVGL EDEEDLLEDL DQALKAHPP SGSHS	
Expression System	E. coli	
Purity	> 95% by SDS-PAGE.	
Conjugation	Unconjugated	
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.	



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX67343-pro Image

3μg Human CTH protein (GTX67343-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).