

Human DPH2 protein, His tag

Cat. No. GTX67359-pro

Species Human

Package

10 µg

Applications

Observed MW (kDa) 53 kDa.

Properties

Form Liquid

Buffer PBS, 20% Glycerol, 1mM DTT

Preservative No preservatives

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration 0.5 mg/ml (Please refer to the vial label for the specific concentration.)

Region/Sequence Full length protein, C-terminal His-Tag; MESMFSSPAE AALQRETGVP GLLTPLDLD GVEYELRVAG FVRDLGCERV ALQFPDQLLG DAVAVAARLE ETTGSKMFIL GDTAYGSCCV DVLGAEQAGA QALIHFGPAC LSPPARPLPV AFVLRQSV A LELCVKAFE A QNPDPKAPV V LLSEPACAH A LEALATLLRP RYLDLLVSSP AFPQPVGSL S PEPMLERFG RRFPLAPGRR LEEYGAFYVG GSKASPD PDL DPDL S RLLLG WAPGQPFSSC CPDTGKTQDE GARAGRLR R RRYLVERARD ARVVGLLAGT LGVAQHREAL AHLRNL TQAA GKRSYVLALG RPTPAKLANF PEVDVFVLLA CPLGALAPQL SGSFFQPILA PCELEAACNP AWPPPGLAPH LTHYADLLPG SPFHVALPPP ESELWETPDV SLITGDLRPP PAWKSSNDHG SLALT PRPQL ELAESSPAAS FLSSRSWQGL EPRLGQTPVT EAVSGRRGIA IAYEDEGSGL EHHHHHHH

Expression System Baculovirus

Purity > 85% by SDS-PAGE.

Endotoxin < 1.0 EU/µg (determined by LAL method)

Conjugation Unconjugated


Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX67359-pro Image

3μg Human DPH2 protein (GTX67359-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).