

Human FABP7 protein

Cat. No. GTX67388-pro

Species Human**Package**
100 µg

Applications

Observed MW (kDa) 15 kDa.

Properties


Form Liquid**Buffer** 25mM Tris-HCl, 10% Glycerol, 1mM DTT, 2mM EDTA**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 1 mg/ml (Please refer to the vial label for the specific concentration.)**Region/Sequence** Full length protein, MVEAFCATWK LTNSQNFDEY MKALGVGFAT RQVGNVTKPT VIISQEGDKV VIRTLSFKN TEISFQLGEE FDETTADDRN CKSVVSLDGD KLVHIQKWDG KETNPFVREIK DGKMVMTLTF GDWAVRHYE KA**Expression System** E. coli**Purity** > 95% by SDS-PAGE.**Conjugation** Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES



GTX67388-pro Image

3µg Human FABP7 protein (GTX67388-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.

For full product information, images and publications, please visit our [website](https://www.genetex.com).

Date 2026 / 10 / 06 Page 1 of 1