

Human FRA2 protein, His tag

Cat. No. GTX67396-pro

Species	Human
---------	-------

Package	50 µg
---------	-------

Applications

Observed MW (kDa)	35 kDa.
-------------------	---------

Properties

Form	Liquid
------	--------

Buffer	20mM Tris-HCl, 10% Glycerol
--------	-----------------------------

Preservative	No preservatives
--------------	------------------

Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
---------	--

Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
---------------	--

Region/Sequence

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMYQDYPG NFDTSSRGSS GSPAHAESYS
SGGGGQKFR VDMPGSGSAF IPTINAITTS QDLQWMVQPT VITMSNPYP RSHPYSLPG LASVPGHMAL PRPGVIKTIG
TTVGRRRRDE QLSPEEEER RIRRRNKLA AAKCRNRRRE LTEKLQAEET EEEEEKSLQ KEIAELQKEK EKLEFMLVAH GPVCKISPEE
RRSPAPGLQ PMRSGGGSVG AVVVKQEPLE EDSPSSSSAG LDKAQRSVIK PISIAGGFYD EEPLHTPIVV TSTPAVTPGT
SNLVFTYPSV LEQESPASPS ESCSKAHRRS SSSGDQSSDS LNSPTLLAL

Expression System	E. coli
-------------------	---------

Purity	> 85% by SDS-PAGE.
--------	--------------------

Conjugation	Unconjugated
-------------	--------------

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX67396-pro Image

3μg Human FRA2 protein (GTX67396-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).