

Human GAD67 protein, His tag

Cat. No. GTX67399-pro

Species Human**Package**
100 µg

Applications

Observed MW (kDa) 67 kDa.

Properties

Form Liquid**Buffer** 20mM Tris-HCl, 10% Glycerol**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 0.5 mg/ml (Please refer to the vial label for the specific concentration.)**Region/Sequence**

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MGSMASSTPS SSATSSNAGA DPNTTNLRPT
TYDTWCGVAH GCTRKLGLKI CGFLQRTNSL EESRLVSAF KERQSSKNLL SCENS DRDAR FRR TETDFSN LFARDLLPAK
NGEEQTVQFL LEVDILLNY VRKTFDRSTK VLD FHHPHQL LEGMEGFNLE LSDHPESLEQ ILVDCRDTLK YGVRTGHPRF
FNQLSTGLDI IGLAGEWLTS TANTNMFTYE IAPVFVLM EQITLKKMREIV GWSSKDGDGI FSPGGAISNM YSIMAARYKY
FPEVKTGMA AVPKLVLFST EQSHYSIKKA GAALGFSTDN VILIKCNERG KIIPADFEAK ILEAKQKGYV PFYVNATAGT TVYGAFDPIQ
EIADICEKYN LWLHVDAAWG GLLMSRKHRL HKLNGIERAN SVTWNPHKMM GVLLQCSAIL VKEKGILQGC NQMCAGYLFQ
PDKQYDVSVD TGDKAIQCGR HVDIFKFWLM WKAKGTVGFE NQINKCLELA EYLYAKIKNR EEFEMVFNGE PEHTNVCFWY
IPQSLRGVPD SPQRREKLHK VAPKIKALMM ESGTTMVGYQ PQGDKANFFR MVISNPAATQ SDIDFLIEEI ERLGQDL

Expression System E. coli**Purity** > 80% by SDS-PAGE.**Conjugation** Unconjugated**Note**


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX67399-pro Image

3μg Human GAD67 protein (GTX67399-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).