

Human HIP2 protein, His tag

Cat. No. GTX67459-pro

Species	Human
---------	-------

Package

50 µg

Applications

Observed MW (kDa)	22 kDa.
-------------------	---------

Properties

Form	Liquid
------	--------

Buffer	20mM Tris-HCl, 50mM NaCl, 10% Glycerol, 1mM DTT
--------	---

Preservative	No preservatives
--------------	------------------

Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
---------	--

Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
---------------	--

Region/Sequence

Full length protein, N-terminal His-Tag; MRGSHHHHHH GMASMTGGQQ MGRDLYDDDD KDRWGSMANI AVQRIKREFK EVLKSEETSK NQIKVDLVDE NfTELrGEIA GPPDTPYEGG RYQLEIKIPE TYPFNPPKVR FITKIWHpNI SSVTGAICLD ILKDQWAAAM TLRTVLLSLQ ALLAAAPDD PQDAVVANQY KQNPEMFkQT ARLWAHVYAG APVSSPEYTK KIENLCAMGF DRNAVIVALs SKSWDVETAT ELLLSN

Expression System	E. coli
-------------------	---------

Purity	> 95% by SDS-PAGE.
--------	--------------------

Conjugation	Unconjugated
-------------	--------------

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX67459-pro Image

3μg Human HIP2 protein (GTX67459-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).