

## Human HNMT protein, His tag

Cat. No. GTX67469-pro

|         |       |
|---------|-------|
| Species | Human |
|---------|-------|

|         |        |
|---------|--------|
| Package | 100 µg |
|---------|--------|

## Applications

|                   |         |
|-------------------|---------|
| Observed MW (kDa) | 33 kDa. |
|-------------------|---------|

## Properties

|      |        |
|------|--------|
| Form | Liquid |
|------|--------|

|        |                             |
|--------|-----------------------------|
| Buffer | 20mM Tris-HCl, 10% Glycerol |
|--------|-----------------------------|

|              |                  |
|--------------|------------------|
| Preservative | No preservatives |
|--------------|------------------|

|         |                                                                                                                                                                                                                            |
|---------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Storage | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
|---------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|               |                                                                          |
|---------------|--------------------------------------------------------------------------|
| Concentration | 1 mg/ml (Please refer to the vial label for the specific concentration.) |
|---------------|--------------------------------------------------------------------------|

## Region/Sequence

Full length protein, N-terminal His-Tag; MRGSHHHHHH GMASMTGGQQ MGRDLYDDDD KDRWGSMAS MRSLFSDHGK YVESFRRFLN HSTEHQCMQE FMDKKLPGL GRIGDTKSEI KLSIGGGAG EIDLQILSKV QAQYPGVCIN NEVVEPSAEQ IAKYKELVAK TSNLENVKFA WHKETSSEYQ SRMLEKKELQ KWDFIHMIQM LYYVKDIPAT LKFFHSLTGT NAKMLIIVVS GSSGWDKLWK KYGSRFPQDD LCQYITSDDL TQMLDNLGLK YECYDLLSTM DISDCFIDGD ENGDLLWDFL TETCNFNATA PPDRLAELGK DLQEPEFSK KEGKVLFNNT LSFIVIA

|                   |         |
|-------------------|---------|
| Expression System | E. coli |
|-------------------|---------|

|        |                    |
|--------|--------------------|
| Purity | > 95% by SDS-PAGE. |
|--------|--------------------|

|             |              |
|-------------|--------------|
| Conjugation | Unconjugated |
|-------------|--------------|

## Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

## DATA IMAGES



### GTX67469-pro Image

3μg Human HNMT protein (GTX67469-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).