

Human HPD protein, His tag

Cat. No. GTX67476-pro

Species Human**Package**

50 µg

Applications

Observed MW (kDa) 45 kDa.

Properties

Form Liquid**Buffer** 20mM Tris-HCl, 50mM NaCl, 20% Glycerol, 1mM DTT**Preservative** No preservatives**Storage** Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.**Concentration** 1 mg/ml (Please refer to the vial label for the specific concentration.)**Region/Sequence**

Full length protein, N-terminal His-Tag; MGSSHHHHHH SSGLVPRGSH MTTYSKDGAK PERGRFLHFH SVTFWVGNAK
QAASFYCSKM GFEPLAYRGL ETGSREVVSH VIKQGKIVFV LSSALNPWNK EMGDHLVKHG DGVKDIAFEV EDCDYIVQKA
RERGAIMRE PWVEQDKFGK VKFAVLQTYG DTTHTLVEKM NYIGQFLPGY EAPAFMDPLL PKLPKCSLEM IDHIVGNQPD
QEMVSASEWY LKNLQFHRFW SVDDTQVHTE YSSLSIVVA NYEESIKMPI NEPAPGKKKS QIQEYVDYNG GAGVQHIALK
TEDIIAIRH LRERGLEFLS VPSTYYQLR EKLKTAKIKV KENIDALEEL KILVDYDEKG YLLQIFTKPV QDRPTLFLEV IQRHNNHGGFG
AGNFNSLFKA FEEEQNLRGN LTNMETNGVV PGM

Expression System E. coli**Purity** > 90% by SDS-PAGE.**Conjugation** Unconjugated**Note**


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX67476-pro Image

3μg Human HPD protein (GTX67476-pro) by SDS-PAGE under reducing condition and visualized by coomassie blue stain.



For full product information, images and publications, please visit our [website](#).